

WesternBright™ Sirius

Most sensitive substrate for chemiluminescent and chemifluorescent Western blotting



WesternBright Sirius is the most sensitive HRP substrate available from Advansta for chemiluminescent Western blotting. With attomolar sensitivity and a long-lasting signal, WesternBright Sirius allows you to detect bands not visualized with other substrates. High signal-to-noise and a large dynamic range make it ideal for quantifying low-intensity bands. And, chemifluorescence emission allows WesternBright Sirius to be imaged by fluorescence imagers in addition to traditional CCD imagers and film.


WesternBright™ Sirius

Most sensitive substrate for chemiluminescent and chemifluorescent Western blotting

WesternBright Sirius allows you to detect low-abundance proteins, with a strong, long-lasting and quantitative signal.

Chemiluminescent and chemifluorescent detection

WesternBright Sirius is the most sensitive chemiluminescent HRP substrate offered by Advansta, specially engineered to produce an extremely strong signal to detect low-abundance proteins and to work with very low antibody concentrations. The substrate also produces a strong chemifluorescent signal that can be detected with fluorescence imagers.



Sensitive chemifluorescent detection

625 pg
313 pg
156 pg
78 pg
39 pg




WesternBright Sirius provides imaging flexibility, with emitted light that can be detected using fluorescent imaging systems in addition to film and CCD based systems. A slot blot containing serial dilutions of an HRP-labeled secondary antibody was detected with WesternBright Sirius, and chemifluorescence imaged using a Typhoon imager (GE Healthcare).

WesternBright™ Sirius

Most sensitive substrate for chemiluminescent and chemifluorescent Western blotting


Image blots hours after substrate incubation

A slot blot containing a serial dilution of an HRP-conjugated antibody was incubated with WesternBright Sirius, and imaged at the times indicated. After 8 hours, the bands are still easily detected with only a 2 minute exposure.



Low background for better band quantification

WesternBright Sirius provides sensitive detection combined with extremely low background. Duplicate Western blots were conducted to analyze samples of plasma, diluted either 1:100 (left lanes) or 1:600 (right lanes) in PBS. The blots were probed with a goat antisera primary antibody, and an HRP-conjugated secondary antibody. After detection with either WesternBright Sirius or SuperSignal® West Femto (Thermo Scientific), the blots were imaged for 8.3 minutes. The signal-to-noise ratio is shown for each band.



WesternBright™ Sirius

Most sensitive substrate for chemiluminescent and chemifluorescent Western blotting


WesternBright Sirius is designed for sensitive detection and low antibody concentrations. With low background and a strong, long-lasting chemiluminescent and chemifluorescent signal, WesternBright Sirius is the best choice for quantifying low-abundance bands.

Features

- Detect low-abundance proteins
- Low background for high signal-to-noise ratio
- Long-lasting signal for flexible imaging
- Image chemiluminescence by CCD or film, image chemifluorescence with laser-based fluorescence imaging systems.



The highest sensitivity



Slot blots containing serial dilutions of an HRP-conjugated antibody were incubated with WesternBright ECL, WesternBright Quantum or WesternBright Sirius and detected with a 2-minute exposure to film.

Ordering Information

Catalog Number	Product	Size
K-12043-C20	WesternBright™ Sirius Western Blotting HRP Substrate Trial kit size (for 200 cm ² membrane)	20 ml
K-12043-D10	WesternBright™ Sirius Western Blotting HRP Substrate (for 1000 cm ² membrane)	100 ml
K-12043-D20	WesternBright™ Sirius Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml

Related Products

L-07001-010	Background Quenching Sheets	10 sheets
L-07013-100	LucentBlue™ X-ray film, 5 x 7" sheets	100 sheets
L-07014-100	LucentBlue™ X-ray film, 8 x 10" sheets	100 sheets
L-08001-010	Low-Fluorescence PVDF Transfer Membrane, 7 x 9 cm	10 sheets
L-08002-010	Nitrocellulose Transfer Membrane 0.45 µm, 8 x 10 cm	10 sheets
L-08003-010	Nitrocellulose Transfer Membrane 0.22 µm, 8 x 10 cm	10 sheets
R-05071-500	Goat-anti-mouse HRP-conjugated secondary antibody	500 µl
R-05072-500	Goat-anti-rabbit HRP-conjugated secondary antibody	500 µl

Copyright © 2011 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ and LucentBlue™ are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

Advansta Corporation
1455 Adams Drive, Ste. 1160 | Menlo Park, CA 94025
Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com
Product information: www.advansta.com/WesternBright_Sirius

