

GRS Genomic DNA Kit – Blood & Cultured Cells

#GK02.0100 (100 preps) | GK02s (trial size, 4 preps)
(FOR RESEARCH ONLY)



SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM BACTERIA

GRAM-NEGATIVE BACTERIA

- 1) Transfer up to 1×10^9 cultured bacterial cells to a 1.5-ml microcentrifuge tube, centrifuge for 1 minute at 14,000g-16,000g and discard the supernatant
- 2) Resuspend the bacterial pellet in 200 μ l of Buffer BC2 and incubate at room temperature for 5 minutes. Proceed with the lysis step of the cultured cell protocol.

GRAM-POSITIVE BACTERIA

Required (and not included):

*Lysozyme Buffer = 20mg/ml Lysozyme; 20mM Tris-HCl pH 8.0; 2mM EDTA; 1% Triton X-100
(prepare fresh immediately prior to use)*

- 1) Transfer up to 1×10^9 cultured bacterial cells to a 1.5-ml microcentrifuge tube, centrifuge for 1 minute at 14,000g-16,000g and discard the supernatant
- 2) Resuspend the bacterial pellet in 200 μ l of Lysozyme Buffer and incubate at room temperature for 10 minutes. During incubation, invert the tube regularly. Proceed with the lysis step of the cultured cell protocol.