

GRS Genomic DNA Kit – Tissue
#GK03.0100 (100 preps) | GK03s (trial size, 4 preps)
(FOR RESEARCH ONLY)



SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM BONE AND TOOTH

- 1) Crush the bone or teeth into small fragments (<5 cm) and grind to a fine powder by either:
 - a) Place fragments in a mortar and pestle and pour in liquid nitrogen. Grind the bone/teeth while immersed in liquid nitrogen and keep pouring liquid nitrogen to ensure the sample remains immersed at all time.
 - b) Place fragments in a metal blender half-filled with liquid nitrogen. Grind the bone/teeth while completely immersed in liquid nitrogen.
- 2) If sample is less than 100mg, transfer into a 1.5-ml microtube and proceed immediately to step 7. If sample is between 100mg and 5g, proceed with step 3.
- 3) Transfer 100mg – 5g to a 50-ml centrifuge tube and add 40 ml of 0.5M EDTA pH 7.5 and incubate with regular agitation (e.g., rotator) at 37°C for 24 hours.
- 4) Centrifuge the sample at 2,000g for 15 minutes and discard the supernatant. Repeat steps 3 and 4 several times for complete decalcification. This process usually requires 3 to 5 days.
- 5) Wash the pellet with 40ml of sterile water. Centrifuge at 2,000g for 15 minutes, discard the supernatant and repeat this washing step three times.
- 6) Transfer up to 50mg into a 1.5ml microtube and proceed with step 7.
- 7) Add 360µl of Buffer BC2 and 40µl of Proteinase K (10mg/ml) and vortex for 30 seconds. Incubate at 60°C until lysis is complete. During incubation, vortex occasionally to disperse the sample.
- 8) Add 400µl of Buffer TC1 and shake vigorously for 10 to 15 seconds.
- 9) Add 400µl of absolute ethanol and shake immediately for 10 to 15 seconds.
- 10) Place a gDNA plus mini column in a 2-ml collection tube and transfer 600 µl of the mixture to the column. Centrifuge at 14,000 to 16,000g for 1 minute. Discard the flow-through, place the column back in the collection tube and transfer all of the remaining sample mixture to the column.
- 11) Centrifuge at 14,000- 16,000g for 1 minute. Discard the collection tube containing the flow-through and place the column into a new 2-ml collection and proceed with the washing step of the normal tissue protocol.