

GRS Genomic DNA Kit – Tissue
#GK03.0100 (100 preps) | GK03s (trial size, 4 preps)
(FOR RESEARCH ONLY)



SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM BUCCAL SWAB

- 1) Add 500µl of Buffer BC2 and 20µl of Proteinase K (10mg/ml) to a 1.5-ml microcentrifuge tube. Place the buccal swab into the tube and incubate at 60°C for 10 minutes.
- 2) Discard the swab and add 500µl of Buffer TC1 to the lysate. Mix immediately by shaking vigorously. Incubate at 60°C for 10 minutes [At this time, preheat the elution buffer to 60°C for the DNA elution step]
- 3) **[optional; when RNA-free DNA is required]** Allow the mixture to cool to room temperature and add 5µl of RNase A (10mg/ml), vortex, and incubate for 5 minutes at room temperature.
- 4) Add 500µl of absolute ethanol to the lysate and immediately mix by shaking vigorously.
- 5) Place the genomic DNA mini spin column in a 2-ml collection tube and transfer 700 µl of the sample mixture to the column. Centrifuge at 14,000g-16,000g for 1 minute.
- 6) Discard the flow-through and repeat step 5 with the remaining lysate.
- 7) Proceed with the wash step of the tissue protocol.