

**GRS Genomic DNA Kit – Tissue**  
#GK03.0100 (100 preps) | GK03s (trial size, 4 preps)  
(FOR RESEARCH ONLY)



## SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM STOOL

### ANIMAL DNA (including parasites) or BACTERIAL DNA

*Required (and not included):*

*Lysozyme Buffer = 20mg/ml Lysozyme; 20mM Tris-HCl pH 8.0; 2mM EDTA; 1% Triton X-100  
(prepare fresh immediately prior to use)*

- 1) Transfer up to 30mg of stool into a 1.5-ml microcentrifuge tube.
- 2) Add 200µl of Lysozyme Buffer and homogenize by vortexing vigorously and/or grinding using the micropestle.
- 3) Incubate at room temperature for 20 minutes. During incubation, invert the tube regularly.
- 4) Proceed with the lysis step of the tissue protocol (increase the temperature to 70°C or 95°C for difficult samples such as Gram-positive bacteria).

### FUNGAL & YEAST DNA

*Required (and not included):*

*Zymolyase or Lyticase*

*Sorbitol Buffer = 1.2M sorbitol; 10mM CaCl<sub>2</sub>; 0.1M Tris-HCl pH 7.5; 35mM β-mercaptoethanol*

- 1) Transfer up to 30mg of stool into a 1.5-ml microcentrifuge tube.
- 2) Add 600µl of Sorbitol Buffer and 200U of Zymolyase (or Lyticase) and homogenize by vortexing vigorously and/or grinding using the micropestle.
- 3) Incubate at 30°C for 30 minutes. During incubation, invert the tube regularly.
- 4) Harvest spheroplasts by centrifugation at 2,000g for 10 minutes. Remove the supernatant and resuspend the pellet in 200µl of Buffer BC2
- 5) Incubate at room temperature for 5 minutes and proceed with the lysis step of the cultured cell protocol.