

GRS Genomic DNA Kit – Card

#GK25.0100 (100 preps) | GK25s (trial size, 4 preps)
(FOR RESEARCH ONLY)



PRIOR TO USE:

Check if Ethanol was added to Wash Buffer 2, water to Proteinase K and Elution Buffer to Carrier RNA. Then, for each sample to be processed, mix 1 µl of RNA Carrier Solution and 200 µl of Buffer C1 in a suitable RNase-free tube. Vortex briefly and keep at hand for step 2. Preheat Elution Buffer to 60°C.

SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM SALIVA

1. Transfer 1-100µl of Urine to a RNase/DNase-free 1.5ml microtube and centrifuge at 6,000g for 2 minutes. Discard the supernatant and add 500 µl of Elution Buffer to the pellet and vortex for 5 seconds.
(For larger volumes: transfer urine to a 15-ml tube and centrifuge at 6,000g for 2 minutes. Discard the supernatant and add 500µl of Elution Buffer. Vortex for 5 seconds and transfer to a 1.5ml microcentrifuge tube)
2. Centrifuge at 6,000g for 2 minutes, discard the supernatant and add 200 µl of Buffer C2 and 20 µl of Proteinase K (10mg/ml), and mix well by vortexing. Lyse by incubation at 60°C for 30 minutes. Invert the tube regularly. *Note: to increase DNA yield from male urine samples, containing sperm cells, add 20µl of 1M DTT (not included) prior to lysis.*
3. Add 200 µl of Buffer C1 (containing RNA Carrier) and mix well by vortexing. Incubate at 60°C for another 20 minutes. Invert the tube regularly.
4. After incubation, centrifuge briefly and transfer the lysate to a new RNase/DNase-free 1.5ml microtube. Add 200µl of ethanol (96%-100%) and mix by shaking vigorously immediately for 10 seconds. In case precipitate appears, break it up by pipetting
5. Place the genomic DNA mini spin column in a 2-ml collection tube and transfer the sample mixture (including any precipitate if present) to the column. Centrifuge at 14,000g-16,000g for 2 minutes.
6. Discard the flow-through and place the spin column in a new 2-ml collection tube and add 400µl of Wash Buffer 1. Centrifuge at 14,000g-16,000g for 30 seconds and discard the flow-through.
7. Place the spin column back in the collection tube. Add 600µl of Wash Buffer 2 (check if ethanol was added) Centrifuge at 14.000g-16.000g for 30 seconds and discard the flow-through.
8. Place the spin column back in the 2-ml collection tube and centrifuge at 14,000g-16,000g for 3 minutes to dry the column membrane completely and to get rid of trace amounts of ethanol. Discard flow-through and collection tube.
9. Remove the spin column carefully and place into a new 1.5-ml microtube.
10. Pipet 50µl Elution Buffer (pre-heated) directly to the center of the spin column without touching the membrane. Incubate at room temperature for 3 minutes.
Notes: Instead of Elution Buffer, DNA can also be diluted with TE or water; pH should be 8.0-8.5. Standard elution volume is 50µl. Concentration can be increased by using less volume (30µl) or alternatively yield can be increased by using more volume (100µl).
11. Centrifuge for at 14,000g-16,000g for 1 minute to elute purified genomic DNA. Discard the spin column and use DNA immediately or store at -20°C.