

## RNase Xterminator Spray

#GB43.500S (500 ml)  
(FOR RESEARCH ONLY)



---

**Product:** Solution for the removal of RNase, other Enzymes, and DNA from laboratory surfaces

**Applications:** Decontamination of work space and materials

**Contents:** 500ml (spray)

**Appearance:** Aqueous, transparent

**Storage:** Store at room temperature

---

### Usage:

RNase Xterminator Spray is a secondary alcohol-based solution containing anionic surfactants for the inactivation and removal of RNases (and other enzymes) as well as DNA from laboratory surfaces. The solution can be used to decontaminate bench tops, labware, pipettors, and other equipment. It is effective on metal, glass and plastic surfaces. Available as a ready-to-use solution, simply spray on the contaminated area and wipe away from the surface using ultrapure water.